thermo scientific

TriPlus RSH Series and TriPlus 100 LS

Multilingual Safety Guide

31709600 Revision E September 2021

ThermoFisher
SCIENTIFIC

© 2021 Thermo Fisher Scientific Inc. All rights reserved.

Thermo Scientific™, TriPlus™ RSH, TriPlus™ RSH SMART, TRACE™ 1300, TRACE™ 1310, TRACE™ GC Ultra, and FOCUS GC are trademarks of Thermo Fisher Scientific Inc., and its subsidiaries.

Published by Thermo Fisher Scientific S.p.A., Strada Rivoltana, km 4, 20090 Rodano - Milan - Italy
Tel: +39 02 950591 - Fax: +39 02 9505276

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Release history:

First edition, released July 2011 "*Original Instructions*"
Second Edition, released January 2012; Third Edition, released October 2013;
Fourth Edition, released November 2016; Fifth Edition, released September 2021.

General Lab Equipment. Not for Clinical, Patient, or Diagnostic Use.

Safety Information and Warnings for Users of the TriPlus RSH Series/TriPlus 100 LS System

This safety guide raises awareness of potential safety issues and general points for consideration for Thermo Fisher Scientific representatives during installation and repair of the TriPlus RSH Series/TriPlus 100 LS or parts of it (following the life cycle principle), as well as for the end user in the lab during the learning phase and in routine work.

It is important to read this section first before operating the TriPlus RSH Series/TriPlus 100 LS.

Contents

- [Product Warranty and Disclaimer](#)
- [General Considerations](#)
- [Environmental Conditions](#)
- [Electrical Hazards](#)
- [Laser Safety Information](#)
- [Other Hazards](#)
- [Working with Toxic or other Harmful Compounds](#)
- [Biological Hazards](#)
- [Maintenance](#)
- [Disposal](#)
- [Regulatory Compliance](#)
- [FCC Compliance Statement](#)
- [RoHS and WEEE Compliance](#)
- [Safety Alerts and Important Information](#)
- [Commonly Used Symbols](#)

Product Warranty and Disclaimer

- Thermo Fisher Scientific reserves the right to make improvements and/or changes to the product(s) described at any time without prior notice.
- Thermo Fisher Scientific makes no warranty of any kind pertaining to this product, including but not limited to implied warranties of merchantability and suitability for a particular purpose.
- Under no circumstances shall Thermo Fisher Scientific be held liable for any coincidental damage or damages arising as a consequence of or from the use of this document.
- Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine readable form without the prior written permission from Thermo Fisher Scientific, except as permitted under copyright laws.
- Thermo Fisher Scientific acknowledges all trade names and trademarks used as the property of their respective owners.

General Considerations

- The TriPlus RSH Series/TriPlus 100 LS User Manual and related documents must be consulted by the user under all circumstances before a unit is put to use.
- Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
- The user should be aware that if the equipment is used in a manner not specified by the manufacturer, the protective and safety features of the equipment may be impaired.
- The repair of instrument failures caused by operation in a manner not specified by the manufacturer is expressly excluded from the standard warranty and service contract coverage.
- When, for technical reasons, it is necessary to work on instrument parts which may involve a potential hazard (moving parts, components under voltage, etc.) the authorized representative of Thermo Fisher Scientific must be contacted. In general, this type of situation arises when access to the parts is only possible using a tool. Should an operator perform a maintenance operation, he/she must have received proper training to carry out that specific task.

Environmental Conditions

The environmental conditions in which the TriPlus RSH Series/TriPlus 100 LS can be operated are listed in the corresponding TriPlus RSH Series/TriPlus 100 LS User Manual; see the section Specifications.

Electrical Hazards

Every analytical instrument has specific hazards, so be sure to read and comply with the following pre-cautions. They will help ensure the safe, long-term use of your TriPlus RSH Series/TriPlus 100 LS.

The installation category (over-voltage category) for this instrument is Level II. The Level II category pertains to equipment receiving its electrical power from the local level, such as an electrical wall outlet.

Connect the TriPlus RSH Series/TriPlus 100 LS only to instruments complying with IEC 61010 safety regulations.

The power line and the connections between the TriPlus RSH Series/TriPlus 100 LS and other instruments used in the configuration setup of the total analytical system must maintain good electrical grounding. Poor grounding represents a danger for the operator and may seriously affect instrument performance.

Do not connect the TriPlus RSH Series/TriPlus 100 LS to power lines that supply devices of a heavy duty nature, such as motors, refrigerators and other devices that can generate electrical disturbances.

Use only fuses of the type and current rating specified. Do not use repaired fuses and do not short-circuit the fuse holder.

The supplied power cord must be inserted into a power outlet with a protective earth contact (ground). When using an extension cord, make sure that the cord also has an earth contact.

If the supplied power cord does not fit the local electrical socket and a replacement or adapter has to be purchased locally, make sure that only a certified power cord is used. Any power cord used should be certified by the appropriate local authorities.

Pay attention not to leave any cable connecting the TriPlus RSH Series/TriPlus 100 LS and the chromatographic system or the power cord close to heated zones, such as the injector or detector heating blocks or the GC hot air vents.

Always replace any cable showing signs of damage with another one provided by the manufacturer.

Safety regulations must be respected.

Do not change the external or internal grounding connections. Tampering with or disconnecting these connections could endanger you and/or damage the TriPlus RSH Series/TriPlus 100 LS.

The instrument is properly grounded in accordance with these regulations when shipped. You do not need to make any changes to the electrical connections or the instrument's chassis to ensure safe operation.

Do not turn the instrument on if you suspect that it has incurred any kind of electrical damage. Instead, disconnect the power cord and contact a Thermo Fisher Scientific representative for a product evaluation. Do not attempt to use the instrument until it has been evaluated. Electrical damage may have occurred if the TriPlus RSH Series/TriPlus 100 LS shows visible signs of damage, exposure to any liquids or has been transported under severe stress.

Damage can also result if the instrument is stored for prolonged periods under unfavorable conditions (e.g. subjected to heat, moisture, etc.).

Ensure that the power supply/controller unit is always placed in a clean and dry position. Avoid any liquid spills in the vicinity.

Before attempting any type of maintenance work, always disconnect the power cord(s) from the power supply(ies) if optional devices are installed. Capacitors inside the instrument may still be charged even if the instrument is turned off.

To avoid damaging electrical parts, do not disconnect an electrical assembly while power is applied to the TriPlus RSH Series/TriPlus 100 LS. Once the power is turned off, wait approximately 30 seconds before you disconnect an assembly.

The instrument includes a number of integrated circuits. These circuits may be damaged if exposed to excessive line voltage fluctuations and/or power surges or electrostatic charges.

Never try to repair or replace any components of the instrument that are not described in the TriPlus RSH Series/TriPlus 100 LS User Manual without the assistance of a Thermo Fisher Scientific representative.

There are no operator-serviceable or replaceable parts inside the power supply(ies) or in the TriPlus RSH Series/TriPlus 100 LS. If a power supply is not functioning, contact a Thermo Fisher Scientific representative.

Laser Safety Information

The power supplies for the TriPlus RSH Series/TriPlus 100 LS, the Temperature Controlled Drawer have the symbols I/O on the label for the power switch to indicate ON/OFF.

If a Temperature Controlled Drawer is installed in combination with a TriPlus RSH Series/TriPlus 100 LS, then a second power supply is active in the complete system. Turning OFF the two power supplies or pulling the two power cords in an emergency will stop the entire TriPlus RSH Series/TriPlus 100 LS.

It is important that the power supply(ies) be in a location where the power ON | OFF switch is accessible and easy to operate, and where it is possible to unplug the AC power cord from the power supply/wall outlet in case of emergency.

Safety Warning for Laser Class 1 Product

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WARNING The installed Laser device is a Class 1 Laser Product. Class 1 Laser devices are not considered to be hazardous when used for their intended purpose. The following statement is required to comply with US and international regulations.

CAUTION Use of controls, adjustments or performance of procedures other than those specified herein may result in hazardous laser light exposures.

The selected Class 1 Laser for the TriPlus RSH Series/TriPlus 100 LS module Barcode Reader complies with the following regulations:

- 21 CFR1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated July 26, 2001
- EN60825-1:2014
- IEC60825-1:2014

The software contains a built-in safety time limit such that the laser scanning mechanism cannot be operated in AIM mode for more than 5 continuous seconds.

Other Hazards

To avoid injury and possible infection through contamination during TriPlus RSH Series/TriPlus 100 LS operation, keep your hands away from the syringe.

Do not operate the TriPlus RSH Series/TriPlus 100 LS without the safety guard. The safety guard must be installed for safe operation.

Do not place any objects inside the area of the safety guard. Keep away from the area around the safety guard during operation of the TriPlus RSH Series/TriPlus 100 LS.

Danger of crushing to fingers and hands. To avoid injury keep your hands away from moving parts during operation. Turn off the power to the TriPlus RSH Series/TriPlus 100 LS if you need to reach inside a mechanically powered system with moving parts.

To avoid injury, observe safe laboratory practice when handling solvents, changing tubing, or operating the TriPlus RSH Series/TriPlus 100 LS. Know the physical and chemical properties of the solvents you use. See the MSDS (Material Safety Data Sheets) from the manufacturer of the solvents being used.

When using the TriPlus RSH Series/TriPlus 100 LS, follow the generally accepted procedures for quality control and method development.

When using the TriPlus RSH Series/TriPlus 100 LS, follow the generally accepted procedures for quality control and method development.

When using the TriPlus RSH Series/TriPlus 100 LS in the field of chromatographic analysis, if a change in the retention of a particular compound, in the resolution between two compounds, or in peak shape is observed, immediately determine the reason(s) for the changes. Do not rely on the separation results until you determine the cause of a change.

Do not operate on the instrument components that form part of the work area of the TriPlus RSH Series/TriPlus 100 LS when this is in motion.

Use caution when working with any polymer tubing under pressure:

- Always wear eye protection when near pressurized polymer tubing.
- Do not use polymer tubing that has been severely stressed or kinked.
- Do not use polymer tubing, in particular no PEEK or Tefzel tubing when using tetrahydrofuran (THF), dimethylsulfoxide (DMSO), chlorinated organic solvents, concentrated mineral acids such as nitric, phosphoric or sulfuric acids, or any related compounds.

Do not use vials without a sealing cap, or microtiter or deepwell plates without a plate seal. Vapor phase from organic solvents can be hazardous and flammable. Acidic vapor phase can cause corrosion to critical mechanical parts.

When sample vials have to undergo heating and agitation, it is important to consider the glass quality. Use high quality glass only. Remember that depending on the application conditions, high pressure can build up in the vial. Whenever a temperature greater than 60 °C is applied, consider the vapor pressure of the solvent used to ensure that no excessive pressure builds up. This is important when using a temperature above 100 °C and especially at the maximum temperature of 200 °C. Be aware that solid materials can also contain volatile compounds such as water (humidity) which could cause build-up of excess vapor pressure.

Do not reuse headspace vials. During the process of washing the vial, micro-cracks can form which will weaken the glass wall and increase the chances of the vial breaking.

In case of a single fault situation where the temperature control of the Agitator fails, there is the potential danger that the device will heat up in an uncontrolled manner until it reaches the cut-off temperature of the overtemperature fuse, in this case, 240 °C.

Based on this single fault scenario, when working with flammable solvents, the user must ensure that the solvent used has a flash point which is 25 °C higher than the maximum potential temperature come in contact with.

When filling-up a standard reservoir or replacing a solvent such as a washing solvent, remove the solvent reservoir bottle from the system to avoid a possible spill over the instrument. Depending on the physical, chemical or hazardous properties of the solvent, use the appropriate protective measures for handling.

Working with Toxic or other Harmful Compounds

Before using dangerous substances (toxic, harmful, etc.) please read the hazard indications and information reported in the Material Safety Data Sheet (MSDS) supplied by the manufacturer referring to the relevant CAS (Chemical Abstract Service) number. The TriPlus RSH Series/TriPlus 100 LS requires the use of several chemical products with different hazard characteristics, which are present in vials and syringes.

Biological Hazards

Before using these substances or replacing the syringe, please read the hazard indications and information reported in the MSDS supplied by the manufacturer referring to the relevant CAS number.

When preparing the samples, please refer to local regulations for the ventilation conditions of the work room.

All waste materials must be collected and eliminated in compliance with the local regulations and directives in the country where the instrument is used.

In laboratories where samples with potential biological hazards are handled, the user must label any equipment or parts thereof which may become contaminated with biohazardous material. The appropriate warning labels are included with the shipment of the instrument. It is the user's responsibility to label the relevant parts of the instrument.

When working with biohazardous materials, it is the user's responsibility to fulfill the following mandatory requirements:

- Instructions on how to safely handle biohazardous material must be provided.
- Operators must be trained and made aware of the potential dangers.
- Personal protective equipment must be provided.
- Instructions must be provided on what to do in case operators are exposed to aerosols or vapors during normal operation (within the intended use of the equipment) or in case of single fault situations such as a broken vial. The protective measures should consider potential contact with the skin, mouth, nose (respiratory organs) and eyes.
- Instructions for decontamination and safe disposal of the relevant parts must be provided.

It is the user's (operator's) responsibility to handle hazardous chemicals or biological compounds, including (but not limited to) bacterial or viral samples and the associated waste, safely and in accordance with international and local regulations.

Maintenance

Any external cleaning or maintenance must be performed with the TriPlus RSH Series/TriPlus 100 LS turned off and the power cord disconnected. Avoid using solvents and spraying on electrical parts. For the removal of potentially dangerous substances (toxic, harmful, etc.) read the hazard indications and information reported in the MSDS (Material Safety Data Sheet) supplied by the manufacturer referring to the relevant CAS (Chemical Abstract Service) number. Use proper protective gloves.

When working with hazardous materials such as radioactive, biologically hazardous material, etc., it is important to train all operators how to respond in case of spills or contamination.

Depending on the class of hazardous material, the appropriate measures have to be taken immediately. Therefore, the chemicals or solvents needed for decontamination have to be on hand.

Any parts of the equipment which can potentially be contaminated, such as the sample vial rack, syringe tool, wash module, etc., should be cleaned regularly. The waste solvent from cleaning and any hardware which needs to be disposed of has to be properly eliminated with all the necessary precautions, abiding by national and international regulations.

When preparing for decontamination, ensure that the solvent or chemical to be used will not damage or react with the surface, dye (color) of the instrument, table or other nearby objects.

If in doubt, please contact your Thermo Fisher Scientific representative to verify the compatibility of the type or composition of solvents with the TriPlus RSH Series/TriPlus 100 LS.

It is the user's (operator's) responsibility to handle hazardous chemicals or biological compounds, including (but not limited to) bacterial or viral samples and the associated wastes, safely and in accordance with international and local regulations.

Disposal

Do not dispose of this equipment or parts thereof unsorted in municipal waste. Follow local municipal waste regulations for proper disposal provisions to reduce the environmental impact of waste electrical and electronic equipment (WEEE).

European Union customers: Call your local customer service representative responsible for the TriPlus RSH Series/TriPlus 100 LS for complimentary equipment pick-up and recycling.

WARNING The customer has to ensure that the TriPlus RSH Series/TriPlus 100 LS has not been contaminated by any hazardous chemical or biological compounds including (but not limited to) bacteria or viruses.

Any part which had direct contact with the analytical sample must be identified and must undergo an appropriate decontamination procedure prior to shipping for disposal.

Potentially dangerous components are: Syringes, Vials and Well Plates. Any critical parts sent for disposal must be handled according to national laws for hazardous compounds.

The customer and the service engineer/technician are fully responsible for enforcing these requirements. Thermo Fisher Scientific will hold the representative and/or customer responsible if these regulations are not observed.

Regulatory Compliance

Thermo Fisher Scientific performs all complete testing and evaluation of its products to ensure full compliance with applicable domestic and international regulations.

Thermo Fisher Scientific declares, under sole responsibility, that the product as originally delivered complies with the requirements of the following applicable European Directives and carries the CE marking accordingly:

- Low Voltage Directive:2014/35/EC
- EMC Directive:2014/30/EC
- Machinery Directive: 2006/42/EC

... and conforms with the following product standards:

Electromagnetic Compatibility (EMC):

- IEC 61326-1:2012, EN 61326-1:2013
- EN 55011:2016+A1:2017, CISPR 11:2015+A1:2016
- FCC part 15, Subpart B, §15.107(a) and §15.109(a)

Safety:

- IEC 61010-1:2010, EN 61010-1:2010,
- IEC 61010-2-010:2014, EN 61010-2-010:2014
- IEC 61010-2-051:2015, EN 61010-2-051:2015
- IEC 61010-2-081:2015, EN 61010-2-081:2015
- UL 61010-1:2012, CSA C22.2 No. 61010-1:2012

Thermo Fisher Scientific hereby declares that our product has been tested and complies with the requirements of FCC Part 15 Subpart B (Class B).

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interferences, and
2. this device must accept any interference received, including interference that may cause undesired operation.

FCC Compliance Statement

RoHS and WEEE Compliance

Supplier Certificate valid for the following product(s) as delivered by Thermo Fisher Scientific:

TriPlus RSH Series/TriPlus 100 LS including the optional modules:

- Referenced Documents

The following European Union directives are referenced in terms of inquiries, specifications, standard(s) or substance list(s):

- Directive 2011/65/EC on Restriction of the use of certain Hazardous Substances, RoHS
- Directive 2012/19/EU on Waste Electrical and Electronic Equipment, WEEE

- Confirmation

We herewith confirm that the above mentioned product(s) does/do not contain any of the reportable substances (substances being banned or to be declared) as listed in the referenced documents in concentrations above the limits as specified therein, respectively – if there are no concentration limits specified – in concentrations being intentionally introduced.

The design and production of the electrical and electronic equipment takes into account the dismantling and recovery of WEEE and its components and materials.

- WEEE Marking

- The products that are subject to the WEEE Directive shipped after July 1st 2006 (date of enforcement of WEEE Directive by the European Union) are compliant with the WEEE marking requirements.
- See WEEE symbol above, section 'Disposal'. This symbol indicates that the product must not be disposed of with household waste.

- Declaration of Compliance

All TriPlus RSH Series/TriPlus 100 LS including the optional modules comply with RoHS and WEEE requirements. There is no substance to be declared.

- Disclaimer

All information in this declaration is given to the best of our present knowledge and beliefs.

The information given does not constitute a warranty within the meaning of warranty law.

Make sure you follow the precautionary notices presented in this guide and in any documentation such as the TriPlus RSH Series/TriPlus 100 LS User Manual. Safety and other special notices include the following:

Notice on the **Proper Use** of the Thermo Fisher Scientific Instruments. In compliance with international regulations: Use of this instrument in a manner not specified by Thermo Fisher Scientific could impair the safety and protection features of the instrument.

WARNING For your safety, and in compliance with international regulations, the physical handling of this Thermo Fisher Scientific instrument requires a team effort to lift and/or move the equipment. This instrument is too heavy and bulky for one person alone to handle safely. Before lifting or moving the TriPlus RSH Series/TriPlus 100 LS, secure all moving parts first to avoid crushing hands or tipping over the unit.


WARNING This is the general safety warning symbol and safety alert word to prevent actions that could cause personal injury. It highlights a general hazard or signifies the presence of another specific safety symbol. This is usually followed by an appropriate safety precautionary message. When you see a safety alert on your instrument or in the documentation, please carefully read the safety instructions before proceeding.

Note Note that highlights information necessary to prevent damage to software, loss of data, invalid test results or instrument damage; might contain information that is critical for optimal performance of the system or other important information about a task.

Safety Alerts and Important Information

Commonly Used Symbols

Symbol	Description
	Caution, or refer to User Manual
	Caution, Risk of Needle-stick Puncture
	Caution, Hot Surface or High Temperature
	Danger of Crushing to Fingers and Hands
	Class 1 Laser Product
	Biological Hazard
	Direct Current
	Alternating Current
	Protective Conductor Terminal, Ground
	Fuse
	Electrical Power ON. Used with Main TriPlus RSH Series/TriPlus 100 LS Power Supply
	Electrical Power OFF. Used with Main TriPlus RSH Series/TriPlus 100 LS Power Supply
	Caution, Risk of Electrical Shock (High Voltage)
	Disposal, Do not dispose in municipal waste. Follow local waste regulations to reduce electrical and electronic waste (WEEE).

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